

Dog B2M/Beta-2-microglobulin ELISA Kit

Catalog No: #EK19256



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dog B2M/Beta-2-microglobulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Beta-2-microglobulin
Target Name	Beta-2-microglobulin
Accession No.	P19341
Uniprot	P19341

Application Details

Assay Type: Sandwich Detect Range: 7.8-500ng/mL Sensitivity: 3.7ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: B2M

Note: This product is for in vitro research use only