

## Pig B2M/Beta-2-microglobulin ELISA Kit

Catalog No: #EK19257

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig B2M/Beta-2-microglobulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Beta-2-microglobulin
Target Name	Beta-2-microglobulin
Other Names	Lactollin
Accession No.	Q07717
Uniprot	Q07717
GeneID	397033;

## Application Details

Assay Type: Sandwich Detect Range: 3.12-200ng/mL Sensitivity: 1.62ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: B2M

Note: This product is for in vitro research use only