

Dog CTSK/Cathepsin K ELISA Kit

Catalog No: #EK19258

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dog CTSK/Cathepsin K ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Cathepsin K
Target Name	Cathepsin K
Accession No.	Q3ZKN1
Uniprot	Q3ZKN1
GeneID	608843;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.18ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CTSK

Note: This product is for in vitro research use only