

Rat Ctsk/Cathepsin K ELISA Kit

Catalog No: #EK19259

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ctsk/Cathepsin K ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cathepsin K
Target Name	Cathepsin K
Accession No.	O35186
Uniprot	O35186
GeneID	29175;

Application Details

Assay Type: Sandwich Detect Range: 6.25-400 pmol/L Sensitivity: 3.125 pmol/L Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ctsk

Note: This product is for in vitro research use only