

Human CCL15/C-C motif chemokine 15 ELISA Kit

Catalog No: #EK19260

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CCL15/C-C motif chemokine 15 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C motif chemokine 15
Target Name	C-C motif chemokine 15
Other Names	Chemokine CC-2,HCC-2,Leukotactin-1,LKN-1,MIP-1 delta,Macrophage inflammatory protein 5,MIP-5,Mrp-2b,NCC-3,Small-inducible cytokine A15,MIP5,NCC3,SCYA15
Accession No.	Q16663
Uniprot	Q16663
GeneID	6359;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.098ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CCL15

Note: This product is for in vitro research use only