

## Human CANX/Calnexin ELISA Kit

Catalog No: #EK19261

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CANX/Calnexin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Calnexin
Target Name	Calnexin
Other Names	IP90, Major histocompatibility complex class I antigen-binding protein p88,p90
Accession No.	P27824
Uniprot	P27824
GeneID	821;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.17ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CANX

Note: This product is for in vitro research use only