

Mouse Canx/Calnexin ELISA Kit

Catalog No: #EK19262

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Canx/Calnexin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calnexin
Target Name	Calnexin
Accession No.	P35564
Uniprot	P35564
GeneID	12330;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.171ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Canx

Note: This product is for in vitro research use only