

Rat Canx/Calnexin ELISA Kit

Catalog No: #EK19263

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Canx/Calnexin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Calnexin
Target Name	Calnexin
Accession No.	P35565
Uniprot	P35565
GeneID	29144;

Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.317ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Canx

Note: This product is for in vitro research use only