

## Human CTSL/Cathepsin L1 ELISA Kit

Catalog No: #EK19265

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CTSL/Cathepsin L1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cathepsin L1
Target Name	Cathepsin L1
Other Names	Cathepsin L, Major excreted protein, MEP, CTSL1
Accession No.	P07711
Uniprot	P07711
GeneID	1514;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20nmol/L Sensitivity: 0.1nmol/L Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CTSL

Note: This product is for in vitro research use only