

Mouse Ctstl/Cathepsin L1 ELISA Kit

Catalog No: #EK19266

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ctstl/Cathepsin L1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cathepsin L1
Target Name	Cathepsin L1
Other Names	Cathepsin L, Major excreted protein, MEP, p39 cysteine proteinase, Ctstl1
Accession No.	P06797
Uniprot	P06797
GeneID	13039;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.3pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ctstl

Note: This product is for in vitro research use only