

Rat Ctst/Cathepsin L1 ELISA Kit

Catalog No: #EK19267

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ctst/Cathepsin L1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cathepsin L1
Target Name	Cathepsin L1
Other Names	Cathepsin L, Cyclic protein 2, CP-2, Major excreted protein, MEP, Ctst1
Accession No.	P07154
Uniprot	P07154
GeneID	25697;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ctst

Note: This product is for in vitro research use only