

Human BIRC6/Baculoviral IAP repeat-containing protein 6 ELISA Kit

Catalog No: #EK19268

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BIRC6/Baculoviral IAP repeat-containing protein 6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Baculoviral IAP repeat-containing protein 6
Target Name	Baculoviral IAP repeat-containing protein 6
Other Names	BIR repeat-containing ubiquitin-conjugating enzyme, BRUCE, RING-type E3 ubiquitin transferase BIRC6, Ubiquitin-conjugating BIR domain enzyme apollon, APOLLON, KIAA1289
Accession No.	Q9NR09
Uniprot	Q9NR09
GeneID	57448;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.056ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BIRC6

Note: This product is for in vitro research use only