

Mouse Col17a1/Collagen alpha-1(XVII) chain ELISA Kit



Catalog No: #EK19269

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Col17a1/Collagen alpha-1(XVII) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen alpha-1(XVII) chain
Target Name	Collagen alpha-1(XVII) chain
Other Names	180 kDa bullous pemphigoid antigen 2,Bullous pemphigoid antigen 2,Bp180,Bpag2
Accession No.	Q07563
Uniprot	Q07563
GeneID	12821;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Col17a1

Note: This product is for in vitro research use only