

Dog APOH/Beta-2-glycoprotein 1 ELISA Kit

Catalog No: #EK19271

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dog APOH/Beta-2-glycoprotein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Beta-2-glycoprotein 1
Target Name	Beta-2-glycoprotein 1
Other Names	Apolipoprotein H,Apo-H,Beta-2-glycoprotein I,B2GPI,Beta(2)GPI
Accession No.	P33703
Uniprot	P33703
GeneID	403945;

Application Details

Assay Type: Sandwich Detect Range:3.12-200ng/mL Sensitivity:1.38ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:APOH

Note: This product is for in vitro research use only