

Mouse Apoh/Beta-2-glycoprotein 1 ELISA Kit

Catalog No: #EK19272

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Apoh/Beta-2-glycoprotein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Beta-2-glycoprotein 1
Target Name	Beta-2-glycoprotein 1
Other Names	APC inhibitor, Activated protein C-binding protein, Apolipoprotein H, Apo-H, Beta-2-glycoprotein 1, B2GPI, B2gp1
Accession No.	Q01339
Uniprot	Q01339
GeneID	11818;

Application Details

Assay Type: Sandwich Detect Range: 12.5-800ng/ml Sensitivity: 6.28ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Apoh

Note: This product is for in vitro research use only