

## Rat Apoh/Beta-2-glycoprotein 1 ELISA Kit

Catalog No: #EK19273

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Apoh/Beta-2-glycoprotein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Beta-2-glycoprotein 1
Target Name	Beta-2-glycoprotein 1
Other Names	Apolipoprotein H,Apo-H,Beta-2-glycoprotein I,B2GPI
Accession No.	P26644
Uniprot	P26644

## Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/ml Sensitivity: 0.79ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Apoh

Note: This product is for in vitro research use only