

# Mouse Cdc42/Cell division control protein 42 homolog ELISA Kit



Catalog No: #EK19275

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cdc42/Cell division control protein 42 homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cell division control protein 42 homolog
Target Name	Cell division control protein 42 homolog
Other Names	G25K GTP-binding protein
Accession No.	P60766
Uniprot	P60766
GeneID	12540;

## Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.91pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cdc42

Note: This product is for in vitro research use only