

Rat Cdc42/Cell division control protein 42 homolog ELISA Kit



Catalog No: #EK19276

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cdc42/Cell division control protein 42 homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cell division control protein 42 homolog
Target Name	Cell division control protein 42 homolog
Accession No.	Q8CFN2
Uniprot	Q8CFN2
GeneID	64465;

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.92pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdc42

Note: This product is for in vitro research use only