

## Rat Anxa3/Annexin A3 ELISA Kit

Catalog No: #EK19278

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Anxa3/Annexin A3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Annexin A3
Target Name	Annexin A3
Other Names	35-alpha calcimedin,Annexin III,Annexin-3,Lipocortin III,Placental anticoagulant protein III,PAP-III,Anx3
Accession No.	P14669
Uniprot	P14669
GeneID	25291;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Anxa3

Note: This product is for in vitro research use only