

Human ANXA1/Annexin A1 ELISA Kit

Catalog No: #EK19279

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ANXA1/Annexin A1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Annexin A1
Target Name	Annexin A1
Other Names	Annexin I, Annexin-1, Calpactin II, Calpactin-2, Chromobindin-9, Lipocortin I, Phospholipase A2 inhibitory protein, p35, ANX1, LPC1
Accession No.	P04083
Uniprot	P04083
GeneID	301;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.1ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ANXA1

Note: This product is for in vitro research use only