

# Human C1QC/Complement C1q subcomponent subunit C ELISA Kit



Catalog No: #EK19280

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Human C1QC/Complement C1q subcomponent subunit C ELISA Kit          |
| Brief Description  | ELISA Kits  |
| Applications       | ELISA   |
| Species Reactivity | Hu  |
| Specificity        | Natural and recombinant Human Complement C1q subcomponent subunit C |
| Target Name        | Complement C1q subcomponent subunit C                               |
| Other Names        | C1QG  |
| Accession No.      | P02747  |
| Uniprot            | P02747  |
| GeneID             | 714;  |

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 33pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: C1QC

Note: This product is for in vitro research use only