

Rat C1qc/Complement C1q subcomponent subunit C ELISA Kit



Catalog No: #EK19282

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat C1qc/Complement C1q subcomponent subunit C ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Complement C1q subcomponent subunit C
Target Name	Complement C1q subcomponent subunit C
Other Names	C1qg
Accession No.	P31722
Uniprot	P31722
GeneID	362634;

Application Details

Assay Type: Sandwich Detect Range: 12.5-800ng/mL Sensitivity: 10ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C1q

Note: This product is for in vitro research use only