

Human CPN2/Carboxypeptidase N subunit 2 ELISA Kit



Catalog No: #EK19283

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CPN2/Carboxypeptidase N subunit 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Carboxypeptidase N subunit 2
Target Name	Carboxypeptidase N subunit 2
Other Names	Carboxypeptidase N 83 kDa chain, Carboxypeptidase N large subunit, Carboxypeptidase N polypeptide 2, Carboxypeptidase N regulatory subunit, ACBP
Accession No.	P22792
Uniprot	P22792
GeneID	1370;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100U/L Sensitivity: 0.78U/L Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CPN2

Note: This product is for in vitro research use only