

Mouse Cpn2/Carboxypeptidase N subunit 2 ELISA Kit

Catalog No: #EK19284

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cpn2/Carboxypeptidase N subunit 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carboxypeptidase N subunit 2
Target Name	Carboxypeptidase N subunit 2
Other Names	Carboxypeptidase N 83 kDa chain, Carboxypeptidase N large subunit, Carboxypeptidase N polypeptide 2, Carboxypeptidase N regulatory subunit
Accession No.	Q9DBB9
Uniprot	Q9DBB9
GeneID	71756;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.69ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cpn2

Note: This product is for in vitro research use only