

Bovine AVP/Vasopressin-neurophysin 2-copeptin ELISA Kit



Catalog No: #EK19285

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine AVP/Vasopressin-neurophysin 2-copeptin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Vasopressin-neurophysin 2-copeptin
Target Name	Vasopressin-neurophysin 2-copeptin
Other Names	AVP-NPII
Accession No.	P01180
Uniprot	P01180
GeneID	280728;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/ml Sensitivity: 21pg/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AVP

Note: This product is for in vitro research use only