

# Mouse Avp/Vasopressin-neurophysin 2-copeptin ELISA Kit



Catalog No: #EK19286

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Mouse Avp/Vasopressin-neurophysin 2-copeptin ELISA Kit           |
| Brief Description  | ELISA Kits   |
| Applications       | ELISA  |
| Species Reactivity | Ms   |
| Specificity        | Natural and recombinant Mouse Vasopressin-neurophysin 2-copeptin |
| Target Name        | Vasopressin-neurophysin 2-copeptin                               |
| Other Names        | AVP-NP11   |
| Accession No.      | P35455   |
| Uniprot            | P35455   |
| GeneID             | 11998;   |

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.77pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Avp

Note: This product is for in vitro research use only