

Pig AVP/Vasopressin-neurophysin 2-copeptin ELISA Kit



Catalog No: #EK19287

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig AVP/Vasopressin-neurophysin 2-copeptin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Vasopressin-neurophysin 2-copeptin
Target Name	Vasopressin-neurophysin 2-copeptin
Other Names	AVP-NPII
Accession No.	P01183
Uniprot	P01183
GeneID	396995;

Application Details

Assay Type: Sandwich Detect Range: 7.8-500pg/mL Sensitivity: 3.7 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AVP

Note: This product is for in vitro research use only