

## Human CPA2/Carboxypeptidase A2 ELISA Kit

Catalog No: #EK19288

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Human CPA2/Carboxypeptidase A2 ELISA Kit          |
| Brief Description  | ELISA Kits  |
| Applications       | ELISA   |
| Species Reactivity | Hu  |
| Specificity        | Natural and recombinant Human Carboxypeptidase A2 |
| Target Name        | Carboxypeptidase A2                               |
| Accession No.      | P48052  |
| Uniprot            | P48052  |
| GeneID             | 1358;   |

## Application Details

Assay Type: Sandwich Detect Range: 1.56-100U/L Sensitivity: 0.78U/L Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CPA2

Note: This product is for in vitro research use only