

## Mouse Cpa2/Carboxypeptidase A2 ELISA Kit

Catalog No: #EK19289

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cpa2/Carboxypeptidase A2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carboxypeptidase A2
Target Name	Carboxypeptidase A2
Accession No.	Q504N0
Uniprot	Q504N0
GeneID	232680;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 40.3pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cpa2

Note: This product is for in vitro research use only