

Rat Crabp2/Cellular retinoic acid-binding protein 2 ELISA Kit

Catalog No: #EK19291

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Crabp2/Cellular retinoic acid-binding protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cellular retinoic acid-binding protein 2
Target Name	Cellular retinoic acid-binding protein 2
Other Names	Cellular retinoic acid-binding protein II, CRABP-II
Accession No.	P51673
Uniprot	P51673
GeneID	29563;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.065ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Crabp2

Note: This product is for in vitro research use only