

# Rat Cpt1a/Carnitine O-palmitoyltransferase 1, liver isoform ELISA Kit

Catalog No: #EK19292

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cpt1a/Carnitine O-palmitoyltransferase 1, liver isoform ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Carnitine O-palmitoyltransferase 1, liver isoform
Target Name	Carnitine O-palmitoyltransferase 1, liver isoform
Other Names	CPT1-L,Carnitine O-palmitoyltransferase I, liver isoform,CPT I,Carnitine palmitoyltransferase 1A,Cpt-1,Cpt1
Accession No.	P32198
Uniprot	P32198
GeneID	25757;

## Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:7.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cpt1a

Note: This product is for in vitro research use only