

Bovine AKT1/RAC-alpha serine/threonine-protein kinase ELISA Kit



Catalog No: #EK19295

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine AKT1/RAC-alpha serine/threonine-protein kinase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine RAC-alpha serine/threonine-protein kinase
Target Name	RAC-alpha serine/threonine-protein kinase
Other Names	Protein kinase B,PKB,Protein kinase B alpha,PKB alpha,RAC-PK-alpha,PKB
Accession No.	Q01314
Uniprot	Q01314
GeneID	280991;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.099ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AKT1

Note: This product is for in vitro research use only