

Rat Akt1/RAC-alpha serine/threonine-protein kinase ELISA Kit



Catalog No: #EK19296

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Akt1/RAC-alpha serine/threonine-protein kinase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat RAC-alpha serine/threonine-protein kinase
Target Name	RAC-alpha serine/threonine-protein kinase
Other Names	Protein kinase B,PKB,Protein kinase B alpha,PKB alpha,RAC-PK-alpha
Accession No.	P47196
Uniprot	P47196
GeneID	24185;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Akt1

Note: This product is for in vitro research use only