

Bovine CTSB/Cathepsin B ELISA Kit

Catalog No: #EK19297

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CTSB/Cathepsin B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cathepsin B
Target Name	Cathepsin B
Other Names	BCSB
Accession No.	P07688
Uniprot	P07688
GeneID	281105;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.09ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CTSB

Note: This product is for in vitro research use only