

Human AHSP/Alpha-hemoglobin-stabilizing protein ELISA Kit



Catalog No: #EK19298

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AHSP/Alpha-hemoglobin-stabilizing protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Alpha-hemoglobin-stabilizing protein
Target Name	Alpha-hemoglobin-stabilizing protein
Other Names	Erythroid differentiation-related factor,Erythroid-associated factor,EDRF,ERAF
Accession No.	Q9NZD4
Uniprot	Q9NZD4
GeneID	51327;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.056ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AHSP

Note: This product is for in vitro research use only