

# Bovine C-Peptide/C-Peptide ELISA Kit

Catalog No: #EK19299



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine C-Peptide/C-Peptide ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine C-Peptide
Target Name	C-Peptide

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.64pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: C-Peptide

Note: This product is for in vitro research use only