

Mouse Cd55/Complement decay-accelerating factor, GPI-anchored ELISA Kit



Catalog No: #EK19301

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd55/Complement decay-accelerating factor, GPI-anchored ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Complement decay-accelerating factor, GPI-anchored
Target Name	Complement decay-accelerating factor, GPI-anchored
Other Names	DAF-GPI,CD55,Cd55a,Daf,Daf1
Accession No.	Q61475
Uniprot	Q61475
GeneID	13136;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.088ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd55

Note: This product is for in vitro research use only