

Bovine BGLAP/Osteocalcin ELISA Kit

Catalog No: #EK19302

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine BGLAP/Osteocalcin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Osteocalcin
Target Name	Osteocalcin
Other Names	Bone Gla protein,BGP,Gamma-carboxyglutamic acid-containing protein
Accession No.	P02820
Uniprot	P02820
GeneID	281646;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.43ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BGLAP

Note: This product is for in vitro research use only