

## Mouse Csta/Cystatin-A ELISA Kit

Catalog No: #EK19305

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Csta/Cystatin-A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cystatin-A
Target Name	Cystatin-A
Other Names	Cystatin-A1, Stefin-A, Csta1
Accession No.	P56567
Uniprot	P56567
GeneID	209294;

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.065ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Csta

Note: This product is for in vitro research use only