

Rat Csta/Cystatin-A ELISA Kit

Catalog No: #EK19306

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Csta/Cystatin-A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cystatin-A
Target Name	Cystatin-A
Other Names	Cystatin-alpha, Epidermal stefin, Epidermal thiol proteinase inhibitor
Accession No.	P01039
Uniprot	P01039
GeneID	288075;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.052ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Csta

Note: This product is for in vitro research use only