

Bovine CKM/Creatine kinase M-type ELISA Kit

Catalog No: #EK19307

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CKM/Creatine kinase M-type ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Creatine kinase M-type
Target Name	Creatine kinase M-type
Other Names	Creatine kinase M chain,M-CK
Accession No.	Q9XSC6
Uniprot	Q9XSC6
GeneID	286822;

Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.38ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CKM

Note: This product is for in vitro research use only