

Human AMD1/S-adenosylmethionine decarboxylase proenzyme ELISA Kit



Catalog No: #EK19309

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AMD1/S-adenosylmethionine decarboxylase proenzyme ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human S-adenosylmethionine decarboxylase proenzyme
Target Name	S-adenosylmethionine decarboxylase proenzyme
Other Names	AdoMetDC,AMD
Accession No.	P17707
Uniprot	P17707
GeneID	262;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.83ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AMD1

Note: This product is for in vitro research use only