

Mouse Amd1/S-adenosylmethionine decarboxylase proenzyme 1 ELISA Kit



Catalog No: #EK19310

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Amd1/S-adenosylmethionine decarboxylase proenzyme 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse S-adenosylmethionine decarboxylase proenzyme 1
Target Name	S-adenosylmethionine decarboxylase proenzyme 1
Other Names	AdoMetDC 1
Accession No.	P0DMN7
Uniprot	P0DMN7
GeneID	11702;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.9pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Amd1

Note: This product is for in vitro research use only