

# Rat Amd1/S-adenosylmethionine decarboxylase proenzyme ELISA Kit

Catalog No: #EK19311

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Amd1/S-adenosylmethionine decarboxylase proenzyme ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat S-adenosylmethionine decarboxylase proenzyme
Target Name	S-adenosylmethionine decarboxylase proenzyme
Other Names	AdoMetDC
Accession No.	P17708
Uniprot	P17708
GeneID	81640;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Amd1

Note: This product is for in vitro research use only