

## Mouse Col15a1/Collagen alpha-1(XV) chain ELISA Kit

Catalog No: #EK19312

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Col15a1/Collagen alpha-1(XV) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen alpha-1(XV) chain
Target Name	Collagen alpha-1(XV) chain
Accession No.	O35206
Uniprot	O35206
GeneID	12819;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.393ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Col15a1

Note: This product is for in vitro research use only