

# Mouse Cartpt/Cocaine- and amphetamine-regulated transcript protein ELISA Kit

Catalog No: #EK19313

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Mouse Cartpt/Cocaine- and amphetamine-regulated transcript protein ELISA Kit        |
| Brief Description  | ELISA Kits  |
| Applications       | ELISA   |
| Species Reactivity | Ms  |
| Specificity        | Natural and recombinant Mouse Cocaine- and amphetamine-regulated transcript protein |
| Target Name        | Cocaine- and amphetamine-regulated transcript protein                               |
| Other Names        | Cart  |
| Accession No.      | P56388  |
| Uniprot            | P56388  |
| GeneID             | 27220;  |

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.089ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cartpt

Note: This product is for in vitro research use only