

## Dog COL1A1/Collagen alpha-1(I) chain ELISA Kit

Catalog No: #EK19314

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dog COL1A1/Collagen alpha-1(I) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Collagen alpha-1(I) chain
Target Name	Collagen alpha-1(I) chain
Other Names	Alpha-1 type I collagen
Accession No.	Q9XSJ7
Uniprot	Q9XSJ7
GeneID	403651;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: COL1A1

Note: This product is for in vitro research use only