

## Mouse Col2a1/Collagen alpha-1(II) chain ELISA Kit

Catalog No: #EK19315

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Col2a1/Collagen alpha-1(II) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen alpha-1(II) chain
Target Name	Collagen alpha-1(II) chain
Other Names	Alpha-1 type II collagen
Accession No.	P28481
Uniprot	P28481
GeneID	12824;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/ml Sensitivity:0.089ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Col2a1

Note: This product is for in vitro research use only