

Rat Aqp1/Aquaporin-1 ELISA Kit

Catalog No: #EK19316

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Aqp1/Aquaporin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aquaporin-1
Target Name	Aquaporin-1
Other Names	AQP-1,Aquaporin-CHIP,Water channel protein for red blood cells and kidney proximal tubule,Chip28
Accession No.	P29975
Uniprot	P29975
GeneID	25240;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/ml Sensitivity:0.28ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aqp1

Note: This product is for in vitro research use only