

Bovine AQP2/Aquaporin-2 ELISA Kit

Catalog No: #EK19317

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine AQP2/Aquaporin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Aquaporin-2
Target Name	Aquaporin-2
Other Names	AQP-2,ADH water channel,Aquaporin-CD,AQP-CD,Collecting duct water channel protein,WCH-CD,Water channel protein for renal collecting duct
Accession No.	P79099
Uniprot	P79099
GeneID	539870;

Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/ml Sensitivity:0.259ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AQP2

Note: This product is for in vitro research use only